

LKexv 3910

MediLine



Total gross capacity	360 litre
Ext. dimensions in mm (w/d/h)	600/615/1840
Int. dimensions in mm (w/d/h)	440/435/1635
Energy consumption in 24 h	0,865 kWh
Climate rating ¹	SN-ST
Cooling system	—
Defrosting	automatic
Temperature range	+3°C to +16°C
Gastronorm	—
Housing material/colour	
Door/lid material	Steel
Liner Material	moulded polystyrol white
Type of control	electronic controller
Temperature display	external digital
Interior light	—
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Shelves	6
Adjustable Shelving	5
Number of drawers	—
Number of baskets	—
Maximum number of baskets	—
Shelf material	glass
Shelf loading	40 kg
Height of adjustable legs	—
Castors/type of castors	—
Door handle	Ergonomic door handle
Foot pedal opener	—
Self-closing door	yes
Insulation	40 - 40 mm
Defrost water drain	—
Refrigerant	R 600a
Airborne acoustical noise emissions	48 dB(A)
Gross/net weight	74,0/68,0 kg
EAN-Nr.	9005382181898

¹ Climate Ratings: **SN** +10°C to +32°C **N** +16°C to +32°C **ST** +16°C to +38°C **T** +16°C to +43°C

Subject to modification. Timestamp: 26.04.2016

LIEBHERR