

11. Technical Report

Power Supply	90-240 VAC 50-60Hz (see identification plate on bottom of centrifuge)
Power requirement	7W , 30mA
Max. Speed	6,600 rpm (Rotor for 1.5/2.0 ml micro tubes)
	6,300 rpm (Rotor for 8-strip PCR tubes)
Max. Centrifugal force	2,000 rcf (Rotor for 1.5/2.0 ml micro tubes)
	2,200 rcf (Rotor for 8-strip PCR tubes)
Ambient temperature	4-40 °C
Max. relative humidity	75%, no condensing moisture
Degree of pollution	N/A
Overvoltage category	N/A
Dimension	Height: 100 mm
	Depth: 161.5 mm
	Width: 160 mm
Weight	0.76 kg(without rotor)
Fuse	250 V: 1 A

12. Warranty Statement

The warranty for the device is 12 months from the date of purchase.

During this period the device will be repaired or exchanged without charge.

Warranty is void on devices with damaged manufacture's number, damaged seal, or damage due to improper use, storage or transport of device or service by unauthorized persons.

CERTIFICATED OF SALE OF EQUIPMENT

Name of company	
Address	
Model Number	
Description	
Serial Number	
Date of Sale	