

Incubator Hood TH 15



Description

The Incubator Hood TH 15 has a temperature range from +5°C above ambient to +50°C and can be combined with all compact shakers KS 15 and Microplate Shakers TiMix 5 or used as independent unit. The height is sufficient for 1000 ml Erlenmeyer flasks.

Basic Equipment

Coated metal housing with 3 windows, the front door opens upwards

Safety cut-off of the shaker if combined with a Bühler Compact Shaker
Standard temperature controller 5°C above ambient to +50°C

Option

[Extended temperature range up to +60°C](#)

Details

Temperature range: 5°C above ambient to +50°C
Temperature accuracy: + / - 1% of the set value
Inner dimensions (w x d x h): 90 x 410 x 310 mm (height sufficient for 1000 ml Erlenmeyer flasks)
Order no.: 6161 000
Accessory: Shelf TH 15

Technical specifications

Temperature sensor: PT 100
Air circulation: 240 m³/h
Electrical supply: 230 V , 50 / 60 Hz (115 V version on request)
Dimensions (w x d x h): 510 x 450 x 455 mm
Weight: 21 kg