

6 Technical Data

Pipetting speed	max. 10 ml/s
Batteries	Type: rechargeable, AAA, NiMH, 1.2 V, 1000 mAh Charging cycles: 500-1000 (when charging as indicated) Charging time: 10 hours Initial charging should be at least 16 hours. Running time: at least 4000 cycles of aspiration and dispensing of 25 ml.
Mains adapter	Input: 100-240 VAC, 50/60 Hz Output: 9 VDC/500 mA
Materials	Housing: ABS Nose piece: POM Pipette mount: Silicone Filter rubber: Silicone
Dimensions (H x W x D)	145 x 130 x 35 mm
Weight	190 g
Ambient conditions	Operation: 5 – 40°C, max. 80% RH Storage: -10 – 50°C, max. 95% RH